



UNIVERSITY HOSPITAL
UNIKLINIKUM
SALZBURG

Depth – dose volume – target volume in IOeRT


Markus Stana, Klarissa Ellmauer, Sara Vockner,
Ivan Messner, Gerd Fastner, Christoph Gaisberger
and Falk Röder

Department of Radiation Therapy and Radiation Oncology
Paracelsus Medical University
Müllner Hauptstrasse 48
5020 Salzburg, Austria


Overview

- Measuring tissue depth
- Calculating the isodose volume
- Calculating the size of the target volume
- What did we learn?


Different methods to measure tissue thickness



Orientation of CBCT slices for depth measurement



Comparison w. manual CBCT measurement




Depth – Master thesis by Klarissa Ellmauer

Assessing the Benefits of 3D Imaging for Intraoperative Electron Radiotherapy (IOERT): A Quantitative Comparison of Geometrical Information in Breast Cancer Treatment


Submitted by: **Klarissa Ellmauer BSc**
at **Master´s Programme MedTech: Functional Imaging,
Conventional and Ion Radiotherapy**
Wiener Neustadt, August 23th 2024

Supervisor: **Markus Stana**

How the program works




Comparison w. manual CBCT measurement




Calculating the size of the isodose volume

- Method 1: body of rotation from measured isodoses in water
- Method 2: volume of a spheroid
- Method 3: isodose volume from 3D voxel data patient model

Method 1: body of rotation from isodoses measured in water



Method 2: volume of a spheroid




$$V = \frac{4}{3}\pi r^2 z$$

r ... radius of tube

z ... half of tissue depth


Method 3: isodose volume from CBCT



Comparison of methods for $V_{90\%}$ calculation

Energy (MeV)	$V_{90\%}$ (ml)		
	Spheroid	measured profiles	Radiance in water
Tube 4cm			
6	2.4	9.0	5.2
9	5.1	19.0	11.0
12	7.2	27.0	23.0
Tube 6cm			
6	9.6	35.0	21.0
9	15.0	55.0	43.0
12	19.8	71.0	62.0
Tube 8cm			
6	17.1	64.0	49.0
9	26.7	109.0	94.0
12	35.2	136.0	123.0

Comparison of calculated $V_{90\%}$ isodose volumes for 110 patients




Tubus:

- Ø 5cm
- Ø 6cm
- Ø 7cm
- Ø 8cm

Target volume best case scenario: spherical tumor



The smallest possible tumor bed volume is given in case of a spherical tumor, as the surface to volume ratio is minimized in a sphere




Size of target volume

$$V_{TV} = \left(\frac{4}{3}\pi \left(\frac{d}{2} + x \right)^3 \right) - \left(\frac{4}{3}\pi \left(\frac{d}{2} \right)^3 \right)$$

d ... diameter of tumor
 x ... irradiation margin



Calculated Target Volumes for different tumor sizes



diameter of tumor d (mm)	tumor volume (ml)	target volume (ml) for different irradiation margins			
		$x = 5 \text{ mm}$	$x = 10 \text{ mm}$	$x = 15 \text{ mm}$	$x = 20 \text{ mm}$
5	0.07	1.70	8.12	22.38	47.65
10	0.52	3.67	13.61	32.99	64.93
15	1.77	6.41	20.68	45.95	85.35
20	4.19	9.95	29.32	61.26	108.91
25	8.18	14.27	39.53	78.93	135.61
30	14.14	19.37	51.31	98.96	165.46
35	22.45	25.26	64.66	121.34	198.44
40	33.51	31.94	79.59	146.08	234.57

$V_{90\%}$ vs. Target Volume for Tubus with 4cm diameter

6 MeV	9 MeV	12 MeV
9 ml	19 ml	27 ml

Diameter of Tumor (mm)	tumor volume (ml)	target volume (ml) and coverage for different irradiation margins			
		x = 5 mm	x = 10 mm	x = 15 mm	x = 20 mm
5	0.07	1.70	8.12	22.38	47.65
10	0.52	3.67	13.61	32.99	64.93
15	1.77	6.41	20.68	45.95	85.35
20	4.19	9.95	29.32	61.26	108.91
25	8.18	14.27	39.53	78.93	135.61
30	14.14	19.37	51.31	98.96	165.46
35	22.45	25.26	64.66	121.34	198.44
40	33.51	31.94	79.59	146.08	234.57

$V_{90\%}$ vs. Target Volume for Tubus with 5cm diameter

6 MeV	9 MeV	12 MeV
21 ml	33 ml	46 ml

Diameter of Tumor (mm)	tumor volume (ml)	target volume (ml) and coverage for different irradiation margins			
		x = 5 mm	x = 10 mm	x = 15 mm	x = 20 mm
5	0.07	1.70	8.12	22.38	47.65
10	0.52	3.67	13.61	32.99	64.93
15	1.77	6.41	20.68	45.95	85.35
20	4.19	9.95	29.32	61.26	108.91
25	8.18	14.27	39.53	78.93	135.61
30	14.14	19.37	51.31	98.96	165.46
35	22.45	25.26	64.66	121.34	198.44
40	33.51	31.94	79.59	146.08	234.57

$V_{90\%}$ vs. Target Volume for Tubus with 6cm diameter

6 MeV	9 MeV	12 MeV
35 ml	55 ml	71 ml

Diameter of Tumor (mm)	tumor volume (ml)	target volume (ml) and coverage for different irradiation margins			
		x = 5 mm	x = 10 mm	x = 15 mm	x = 20 mm
5	0.07	1.70	8.12	22.38	47.65
10	0.52	3.67	13.61	32.99	64.93
15	1.77	6.41	20.68	45.95	85.35
20	4.19	9.95	29.32	61.26	108.91
25	8.18	14.27	39.53	78.93	135.61
30	14.14	19.37	51.31	98.96	165.46
35	22.45	25.26	64.66	121.34	198.44
40	33.51	31.94	79.59	146.08	234.57

$V_{90\%}$ vs. Target Volume for Tubus with 7cm diameter

6 MeV	9 MeV	12 MeV
48 ml	81 ml	104 ml

Diameter of Tumor (mm)	tumor volume (ml)	target volume (ml) and coverage for different irradiation margins			
		x = 5 mm	x = 10 mm	x = 15 mm	x = 20 mm
5	0.07	1.70	8.12	22.38	47.65
10	0.52	3.67	13.61	32.99	64.93
15	1.77	6.41	20.68	45.95	85.35
20	4.19	9.95	29.32	61.26	108.91
25	8.18	14.27	39.53	78.93	135.61
30	14.14	19.37	51.31	98.96	165.46
35	22.45	25.26	64.66	121.34	198.44
40	33.51	31.94	79.59	146.08	234.57

$V_{90\%}$ vs. Target Volume for Tubus with 8cm diameter

6 MeV	9 MeV	12 MeV
64 ml	109 ml	136 ml

Diameter of Tumor (mm)	tumor volume (ml)	target volume (ml) and coverage for different irradiation margins			
		x = 5 mm	x = 10 mm	x = 15 mm	x = 20 mm
5	0.07	1.70	8.12	22.38	47.65
10	0.52	3.67	13.61	32.99	64.93
15	1.77	6.41	20.68	45.95	85.35
20	4.19	9.95	29.32	61.26	108.91
25	8.18	14.27	39.53	78.93	135.61
30	14.14	19.37	51.31	98.96	165.46
35	22.45	25.26	64.66	121.34	198.44
40	33.51	31.94	79.59	146.08	234.57


Size of target volume considering resection margin

$$V_{TV} = \left(\frac{4}{3}\pi \left(\frac{d}{2} + r + x \right)^3 \right) - \left(\frac{4}{3}\pi \left(\frac{d}{2} + r \right)^3 \right)$$

d ... diameter of tumor

r ... resection margin

x ... irradiation margin



Target Volumes for different tumor sizes and $r = 5\text{mm}$ resection margin

Diameter of Tumor d (mm)	tumor volume (ml)	target volume (ml) for different irradiation margins		
		$x = 5 \text{ mm}$	$x = 10 \text{ mm}$	$x = 15 \text{ mm}$
5	0.07	6.41	20.68	45.95
10	0.52	9.95	29.32	61.26
15	1.77	14.27	39.53	78.93
20	4.19	19.37	51.31	98.96
25	8.18	25.26	64.66	121.34
30	14.14	31.94	79.59	146.08
35	22.45	39.40	96.08	173.18
40	33.51	47.65	114.14	202.63

$V_{90\%}$ vs. TV for 5mm resection margin and a 4cm Tubus

6 MeV	9 MeV	12 MeV
9 ml	19 ml	27 ml

Diameter of Tumor d (mm)	tumor volume (ml)	target volume (ml) and coverage for different irradiation margins		
		$x = 5$ mm	$x = 10$ mm	$x = 15$ mm
5	0.07	6.41	20.68	45.95
10	0.52	9.95	29.32	61.26
15	1.77	14.27	39.53	78.93
20	4.19	19.37	51.31	98.96
25	8.18	25.26	64.66	121.34
30	14.14	31.94	79.59	146.08
35	22.45	39.40	96.08	173.18
40	33.51	47.65	114.14	202.63

$V_{90\%}$ vs. TV for 5mm resection margin and a 5cm Tubus

6 MeV	9 MeV	12 MeV
21 ml	33 ml	46 ml

Diameter of Tumor d (mm)	tumor volume (ml)	target volume (ml) and coverage for different irradiation margins		
		$x = 5$ mm	$x = 10$ mm	$x = 15$ mm
5	0.07	6.41	20.68	45.95
10	0.52	9.95	29.32	61.26
15	1.77	14.27	39.53	78.93
20	4.19	19.37	51.31	98.96
25	8.18	25.26	64.66	121.34
30	14.14	31.94	79.59	146.08
35	22.45	39.40	96.08	173.18
40	33.51	47.65	114.14	202.63

$V_{90\%}$ vs. TV for 5mm resection margin and a 6cm Tubus

6 MeV	9 MeV	12 MeV
35 ml	55 ml	71 ml

Diameter of Tumor d (mm)	tumor volume (ml)	target volume (ml) and coverage for different irradiation margins		
		$x = 5$ mm	$x = 10$ mm	$x = 15$ mm
5	0.07	6.41	20.68	45.95
10	0.52	9.95	29.32	61.26
15	1.77	14.27	39.53	78.93
20	4.19	19.37	51.31	98.96
25	8.18	25.26	64.66	121.34
30	14.14	31.94	79.59	146.08
35	22.45	39.40	96.08	173.18
40	33.51	47.65	114.14	202.63

$V_{90\%}$ vs. TV for 5mm resection margin and a 7cm Tubus

6 MeV	9 MeV	12 MeV
48 ml	81 ml	104 ml

Diameter of Tumor d (mm)	tumor volume (ml)	target volume (ml) and coverage for different irradiation margins		
		$x = 5$ mm	$x = 10$ mm	$x = 15$ mm
5	0.07	6.41	20.68	45.95
10	0.52	9.95	29.32	61.26
15	1.77	14.27	39.53	78.93
20	4.19	19.37	51.31	98.96
25	8.18	25.26	64.66	121.34
30	14.14	31.94	79.59	146.08
35	22.45	39.40	96.08	173.18
40	33.51	47.65	114.14	202.63

$V_{90\%}$ vs. TV for 5mm resection margin and a 8cm Tubus

6 MeV	9 MeV	12 MeV
64 ml	109 ml	136 ml

Diameter of Tumor d (mm)	tumor volume (ml)	target volume (ml) and coverage for different irradiation margins		
		$x = 5$ mm	$x = 10$ mm	$x = 15$ mm
5	0.07	6.41	20.68	45.95
10	0.52	9.95	29.32	61.26
15	1.77	14.27	39.53	78.93
20	4.19	19.37	51.31	98.96
25	8.18	25.26	64.66	121.34
30	14.14	31.94	79.59	146.08
35	22.45	39.40	96.08	173.18
40	33.51	47.65	114.14	202.63

Target Volumes for different tumor sizes and $r = 10\text{mm}$ resection margin

Diameter of Tumor d (mm)	tumor volume (ml)	target volume (ml) for different irradiation margins	
		$x = 5 \text{ mm}$	$x = 10 \text{ mm}$
5	0.07	14.27	39.53
10	0.52	19.37	51.31
15	1.77	25.26	64.66
20	4.19	31.94	79.59
25	8.18	39.40	96.08
30	14.14	47.65	114.14
35	22.45	56.68	133.78
40	33.51	66.50	154.99

$V_{90\%}$ vs. TV for 10mm resection margin and a 4cm Tubus

6 MeV	9 MeV	12 MeV
9 ml	19 ml	27 ml

Diameter of Tumor d (mm)	tumor volume (ml)	target volume (ml) and coverage for different irradiation margins			
		$x = 5 \text{ mm}$		$x = 10 \text{ mm}$	
5	0.07			14.27	39.53
10	0.52			19.37	51.31
15	1.77			25.26	64.66
20	4.19			31.94	79.59
25	8.18			39.40	96.08
30	14.14			47.65	114.14
35	22.45			56.68	133.78
40	33.51			66.50	154.99

$V_{90\%}$ vs. TV for 10mm resection margin and a 5cm Tubus

6 MeV	9 MeV	12 MeV
21 ml	33 ml	46 ml

Diameter of Tumor d (mm)	tumor volume (ml)	target volume (ml) and coverage for different irradiation margins		
		x = 5 mm	x = 10 mm	x = 15 mm
5	0.07	14.27	39.53	
10	0.52	19.37	51.31	
15	1.77	25.26	64.66	
20	4.19	31.94	79.59	
25	8.18	39.40	96.08	
30	14.14	47.65	114.14	
35	22.45	56.68	133.78	
40	33.51	66.50	154.99	

$V_{90\%}$ vs. TV for 10mm resection margin and a 6cm Tubus

6 MeV	9 MeV	12 MeV
35 ml	55 ml	71 ml

Diameter of Tumor d (mm)	tumor volume (ml)	target volume (ml) and coverage for different irradiation margins		
		x = 5 mm	x = 10 mm	x = 15 mm
5	0.07	14.27	39.53	
10	0.52	19.37	51.31	
15	1.77	25.26	64.66	
20	4.19	31.94	79.59	
25	8.18	39.40	96.08	
30	14.14	47.65	114.14	
35	22.45	56.68	133.78	
40	33.51	66.50	154.99	

$V_{90\%}$ vs. TV for 10mm resection margin and a 7cm Tubus

6 MeV	9 MeV	12 MeV
48 ml	81 ml	104 ml

Diameter of Tumor d (mm)	tumor volume (ml)	target volume (ml) and coverage for different irradiation margins		
		x = 5 mm	x = 10 mm	x = 15 mm
5	0.07	14.27	39.53	
10	0.52	19.37	51.31	
15	1.77	25.26	64.66	
20	4.19	31.94	79.59	
25	8.18	39.40	96.08	
30	14.14	47.65	114.14	
35	22.45	56.68	133.78	
40	33.51	66.50	154.99	

$V_{90\%}$ vs. TV for 10mm resection margin and a 8cm Tubus

6 MeV	9 MeV	12 MeV
64 ml	109 ml	136 ml

Diameter of Tumor d (mm)	tumor volume (ml)	target volume (ml) and coverage for different irradiation margins		
			$x = 5$ mm	$x = 10$ mm
5	0.07		14.27	39.53
10	0.52		19.37	51.31
15	1.77		25.26	64.66
20	4.19		31.94	79.59
25	8.18		39.40	96.08
30	14.14		47.65	114.14
35	22.45		56.68	133.78
40	33.51		66.50	154.99

What we (I) have learned

- Maximizing the size of the irradiated volume is probably a good idea
- Therefore the distance to the nearest organ at risk has to be measured exactly
- Knowing the size of the target volume and comparing it with the irradiated volume (a certain isodose volume) should be considered as a plan optimization parameter